

Otoferlin Polyclonal Antibody FITC Conjugated

Catalog No: #C91271FITC



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Otoferlin Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1601-1700/1997
Conjugates	FITC
Target Name	Otoferlin
Other Names	Otoferlin
Accession No.	Uniprot:Q9HC10Gene ID:9381
Uniprot	Q9HC10
GeneID	9381
Calculated MW	227
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only