

METTL3 Polyclonal Antibody APC Conjugated

Catalog No: #C91276APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	METTL3 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	2-100/580
Conjugates	APC
Target Name	METTL3
Other Names	METTL3
Accession No.	Uniprot:Q86U44Gene ID:56339
Uniprot	Q86U44
GeneID	56339
Calculated MW	64
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only