

TET3 Polyclonal Antibody FITC Conjugated

Catalog No: #C91278FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TET3 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1561-1660/1660
Conjugates	FITC
Target Name	TET3
Other Names	TET3
Accession No.	Uniprot:O43151Gene ID:200424
Uniprot	O43151
GeneID	200424
Calculated MW	179
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only