

## FDC-SP Polyclonal Antibody Biotin Conjugated

Catalog No: #C91281Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FDC-SP Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-P IHC-F IF
Species Reactivity	Hu Ms Rat
Immunogen Type	18-85/85
Conjugates	Bio
Target Name	FDC-SP
Other Names	FDC-SP
Accession No.	Uniprot:Q5R1M6Gene ID:260436
Uniprot	Q5R1M6
GeneID	260436
Calculated MW	7.8
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only