

phospho-MRCL3 (Thr18 / Ser19) Polyclonal Antibody HRP Conjugated



Catalog No: #C91287HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-MRCL3 (Thr18 / Ser19) Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-P IHC-F ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	RA(p-T)(p-S)N
Conjugates	HRP
Target Name	phospho-MRCL3 (Thr18 / Ser19)
Other Names	phospho-MRCL3 (Thr18 / Ser19)
Accession No.	Uniprot:P19105Gene ID:10627
Uniprot	P19105
GeneID	10627
Calculated MW	20
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only