

CD71/Transferrin receptor Polyclonal Antibody Biotin Conjugated



Catalog No: #C91291Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---------------------------------------------------------------------------------|
| Product Name | CD71/Transferrin receptor Polyclonal Antibody Biotin Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | WB |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 21-120/760 |
| Conjugates | Bio |
| Target Name | CD71/Transferrin receptor |
| Other Names | CD71/Transferrin receptor |
| Accession No. | Uniprot:P02786Gene ID:7037 |
| Uniprot | P02786 |
| GeneID | 7037 |
| Calculated MW | 74/85 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB=1:50-200

Note: This product is for in vitro research use only