

CK14/17/42/10 Polyclonal Antibody FITC Conjugated

Catalog No: #C91297FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CK14/17/42/10 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	251-350/472
Conjugates	FITC
Target Name	CK14/17/42/10
Other Names	CK14/17/42/10
Accession No.	Uniprot:P02533Gene ID:3861
Uniprot	P02533
GeneID	3861
Calculated MW	52
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only