

LDHA Polyclonal Antibody FITC Conjugated

Catalog No: #C91307FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LDHA Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	261-361/361
Conjugates	FITC
Target Name	LDHA
Other Names	LDHA
Accession No.	Uniprot:P00338Gene ID:3939
Uniprot	P00338
GeneID	3939
Calculated MW	37
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only