

## PPARGC1A Polyclonal Antibody FITC Conjugated

Catalog No: #C91321FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PPARGC1A Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	601-700/798
Conjugates	FITC
Target Name	PPARGC1A
Other Names	PPARGC1A
Accession No.	Uniprot:Q9UBK2Gene ID:10891
Uniprot	Q9UBK2
GeneID	10891
Calculated MW	88
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only