

Lamin B Polyclonal Antibody PE Conjugated

Catalog No: #C91327PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lamin B Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	361-460/586
Conjugates	PE
Target Name	Lamin B
Other Names	Lamin B
Accession No.	Uniprot:P20700Gene ID:4001
Uniprot	P20700
GeneID	4001
Calculated MW	64
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only