

PD-1 Polyclonal Antibody PE-Cy7 Conjugated

Catalog No: #C91341PE-Cy7

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PD-1 Polyclonal Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-288/288
Conjugates	PE-Cy7
Target Name	PD-1
Other Names	PD-1
Accession No.	Uniprot:Q15116Gene ID:5133
Uniprot	Q15116
GeneID	5133
Calculated MW	32
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only