

Matrilin 2 Polyclonal Antibody FITC Conjugated

Catalog No: #C91342FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Matrilin 2 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	601-700/956
Conjugates	FITC
Target Name	Matrilin 2
Other Names	Matrilin 2
Accession No.	Uniprot:O00339Gene ID:4147
Uniprot	O00339
GeneID	4147
Calculated MW	104
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:200-800

Note: This product is for in vitro research use only