

# Serotonin transporter Polyclonal Antibody Biotin Conjugated

Catalog No: #C91352Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serotonin transporter Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P
Species Reactivity	Hu Ms Rat
Immunogen Type	401-500/630
Conjugates	Bio
Target Name	Serotonin transporter
Other Names	Serotonin transporter
Accession No.	Uniprot:P31645Gene ID:6532
Uniprot	P31645
GeneID	6532
Calculated MW	70
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-P=1:50-200

Note: This product is for in vitro research use only