

SOCS2 Polyclonal Antibody FITC Conjugated

Catalog No: #C91354FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SOCS2 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1-100/198
Conjugates	FITC
Target Name	SOCS2
Other Names	SOCS2
Accession No.	Uniprot:O14508Gene ID:8835
Uniprot	O14508
GeneID	8835
Calculated MW	23
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only