

## SSTR4 Polyclonal Antibody FITC Conjugated

Catalog No: #C91355FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	SSTR4 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/388
Conjugates	FITC
Target Name	SSTR4
Other Names	SSTR4
Accession No.	Uniprot:P31391Gene ID:6754
Uniprot	P31391
GeneID	6754
Calculated MW	42
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only