

Syndecan 2 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91358Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Syndecan 2 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-F IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-201/201
Conjugates	Bio
Target Name	Syndecan 2
Other Names	Syndecan 2
Accession No.	Uniprot:P34741Gene ID:6383
Uniprot	P34741
GeneID	6383
Calculated MW	20
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only