

TGF beta Receptor III Polyclonal Antibody Biotin Conjugated

Catalog No: #C91359Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TGF beta Receptor III Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F IF
Species Reactivity	Hu Ms Rat
Immunogen Type	701-850/850
Conjugates	Bio
Target Name	TGF beta Receptor III
Other Names	TGF beta Receptor III
Accession No.	Uniprot:Q03167Gene ID:7049
Uniprot	Q03167
GeneID	7049
Calculated MW	94
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only