

# TGF beta Receptor III Polyclonal Antibody PE-Cy7 Conjugated

Catalog No: #C91359PE-Cy7

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TGF beta Receptor III Polyclonal Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	701-850/850
Conjugates	PE-Cy7
Target Name	TGF beta Receptor III
Other Names	TGF beta Receptor III
Accession No.	Uniprot:Q03167Gene ID:7049
Uniprot	Q03167
GeneID	7049
Calculated MW	94
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only