

CHRM5 Polyclonal Antibody FITC Conjugated

Catalog No: #C91360FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CHRM5 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	461-532/532
Conjugates	FITC
Target Name	CHRM5
Other Names	CHRM5
Accession No.	Uniprot:P08912Gene ID:1133
Uniprot	P08912
GeneID	1133
Calculated MW	60
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only