

phospho-Nephrin (Tyr1217) Polyclonal Antibody Biotin Conjugated



Catalog No: #C91366Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | phospho-Nephrin (Tyr1217) Polyclonal Antibody Biotin Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | WB IHC-P IHC-F ICC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | GI(p-Y)DQ |
| Conjugates | Bio |
| Target Name | phospho-Nephrin (Tyr1217) |
| Other Names | phospho-Nephrin (Tyr1217) |
| Accession No. | Uniprot:O60500Gene ID:4868 |
| Uniprot | O60500 |
| GeneID | 4868 |
| Calculated MW | 136 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only