

Nesprin 3 Polyclonal Antibody FITC Conjugated

Catalog No: #C91369FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Nesprin 3 Polyclonal Antibody FITC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | ICC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 201-300/975 |
| Conjugates | FITC |
| Target Name | Nesprin 3 |
| Other Names | Nesprin 3 |
| Accession No. | Uniprot:Q6ZMZ3Gene ID:161176 |
| Uniprot | Q6ZMZ3 |
| GeneID | 161176 |
| Calculated MW | 112 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only