

## NISCH Polyclonal Antibody FITC Conjugated

Catalog No: #C91372FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NISCH Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Ms Rat
Immunogen Type	1251-1350/1504
Conjugates	FITC
Target Name	NISCH
Other Names	NISCH
Accession No.	Uniprot:Q9Y2I1 Gene ID:11188
Uniprot	Q9Y2I1
GeneID	11188
Calculated MW	167
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only