

## uPA Receptor Polyclonal Antibody Biotin Conjugated

Catalog No: #C91377Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	uPA Receptor Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/335
Conjugates	Bio
Target Name	uPA Receptor
Other Names	uPA Receptor
Accession No.	Uniprot:Q03405Gene ID:5329
Uniprot	Q03405
GeneID	5329
Calculated MW	31
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only