

# Norovirus Capsid protein VP1 Polyclonal Antibody HRP Conjugated



Catalog No: #C91382HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Norovirus Capsid protein VP1 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC
Species Reactivity	Norovirus.
Immunogen Type	341-440/530
Conjugates	HRP
Target Name	Norovirus Capsid protein VP1
Other Names	Norovirus Capsid protein VP1
Accession No.	Uniprot:Q83884Gene ID:1491972
Uniprot	Q83884
GeneID	1491972
Calculated MW	58
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only