

SAA1/Serum Amyloid A Polyclonal Antibody FITC Conjugated



Catalog No: #C91386FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SAA1/Serum Amyloid A Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu
Immunogen Type	41-94/122
Conjugates	FITC
Target Name	SAA1/Serum Amyloid A
Other Names	SAA1/Serum Amyloid A
Accession No.	Uniprot:P0DJ18Gene ID:6288
Uniprot	P0DJ18
GeneID	6288
Calculated MW	14
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only