

## Selenoprotein O Polyclonal Antibody HRP Conjugated

Catalog No: #C91397HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Selenoprotein O Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC
Species Reactivity	Hu Ms
Immunogen Type	301-400/669
Conjugates	HRP
Target Name	Selenoprotein O
Other Names	Selenoprotein O
Accession No.	Uniprot:Q9BVL4Gene ID:83642
Uniprot	Q9BVL4
GeneID	83642
Calculated MW	73
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only