

Lipoprotein lipase Polyclonal Antibody Biotin Conjugated

Catalog No: #C91404Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lipoprotein lipase Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F IF
Species Reactivity	Hu Ms Rat
Immunogen Type	401-475/475
Conjugates	Bio
Target Name	Lipoprotein lipase
Other Names	Lipoprotein lipase
Accession No.	Uniprot:P06858Gene ID:4023
Uniprot	P06858
GeneID	4023
Calculated MW	52
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only