

SH3GL2/Endophilin A1 Polyclonal Antibody Bio Conjugated

Catalog No: #C91406Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SH3GL2/Endophilin A1 Polyclonal Antibody Bio Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	Bio
Target Name	SH3GL2/Endophilin A1
Other Names	SH3GL2/Endophilin A1
Accession No.	Uniprot:Q99962Gene ID:512450
Uniprot	Q99962
GeneID	512450
Calculated MW	40
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only