

phospho-NFAT2 (Ser233) Polyclonal Antibody HRP Conjugated



Catalog No: #C91407HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-NFAT2 (Ser233) Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	LG(p-S)SP
Conjugates	HRP
Target Name	phospho-NFAT2 (Ser233)
Other Names	phospho-NFAT2 (Ser233)
Accession No.	Uniprot:O95644 Gene ID:4772
Uniprot	O95644
GeneID	4772
Calculated MW	101
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only