

SLC34A1 Polyclonal Antibody FITC Conjugated

Catalog No: #C91409FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLC34A1 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/639
Conjugates	FITC
Target Name	SLC34A1
Other Names	SLC34A1
Accession No.	Uniprot:Q06495Gene ID:6569
Uniprot	Q06495
GeneID	6569
Calculated MW	69
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only