

SLC25A38 Polyclonal Antibody FITC Conjugated

Catalog No: #C91410FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLC25A38 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/304
Conjugates	FITC
Target Name	SLC25A38
Other Names	SLC25A38
Accession No.	Uniprot:Q96DW6Gene ID:54977
Uniprot	Q96DW6
GeneID	54977
Calculated MW	34
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:200-800

Note: This product is for in vitro research use only