

## RUNX2 Polyclonal Antibody HRP Conjugated

Catalog No: #C91418HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RUNX2 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/521
Conjugates	HRP
Target Name	RUNX2
Other Names	RUNX2
Accession No.	Uniprot:Q13950Gene ID:860
Uniprot	Q13950
GeneID	860
Calculated MW	57(hu)/67(mo, rat)
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only