

RUNX2 Polyclonal Antibody FITC Conjugated

Catalog No: #C91418FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RUNX2 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/521
Conjugates	FITC
Target Name	RUNX2
Other Names	RUNX2
Accession No.	Uniprot:Q13950Gene ID:860
Uniprot	Q13950
GeneID	860
Calculated MW	57(hu)/67(mo, rat)
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only