

# Placental alkaline phosphatase (PLAP) Polyclonal Antibody PE Conjugated

Catalog No: #C91420PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Placental alkaline phosphatase (PLAP) Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms
Immunogen Type	41-140/535
Conjugates	PE
Target Name	Placental alkaline phosphatase (PLAP)
Other Names	Placental alkaline phosphatase (PLAP)
Accession No.	Uniprot:P05187Gene ID:250
Uniprot	P05187
GeneID	250
Calculated MW	53
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only