

PCV2 Cap protein Polyclonal Antibody HRP Conjugated

Catalog No: #C91423HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PCV2 Cap protein Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	PCV2
Immunogen Type	full length protein
Conjugates	HRP
Target Name	PCV2 Cap protein
Other Names	PCV2 Cap protein
Accession No.	Uniprot:O56129Gene ID:2943185
Uniprot	O56129
GeneID	2943185
Calculated MW	28
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000

Note: This product is for in vitro research use only