

# H1N1 Hemagglutinin 1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91425Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	H1N1 Hemagglutinin 1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F IF
Species Reactivity	Influenza A virus H1N1
Immunogen Type	151-250/566
Conjugates	Bio
Target Name	H1N1 Hemagglutinin 1
Other Names	H1N1 Hemagglutinin 1
Accession No.	Uniprot:P03452Gene ID:23308115
Uniprot	P03452
GeneID	23308115
Calculated MW	63
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only