

CSFV E2 protein Polyclonal Antibody bs-20116G-Biotin Conjugated

Catalog No: #C91432bs-20116G-Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CSFV E2 protein Polyclonal Antibody bs-20116G-Biotin Conjugated
Purification	affinity purified by Protein A
Applications	WB IF
Species Reactivity	CSFV
Conjugates	bs-20116G-Bio
Target Name	CSFV E2 protein
Other Names	CSFV E2 protein
Accession No.	Uniprot:P19712Gene ID:920146
Uniprot	P19712
GeneID	920146
Calculated MW	42
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only