

SCN1A Polyclonal Antibody Biotin Conjugated

Catalog No: #C91433Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SCN1A Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	601-700/2009
Conjugates	Bio
Target Name	SCN1A
Other Names	SCN1A
Accession No.	Uniprot:P35498Gene ID:6323
Uniprot	P35498
GeneID	6323
Calculated MW	229
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only