

# phospho-Nrf2 (Ser40) Polyclonal Antibody Biotin Conjugated

Catalog No: #C91434Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	phospho-Nrf2 (Ser40) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F IF
Species Reactivity	Hu Ms Rat
Immunogen Type	DF(p-S)QR
Conjugates	Bio
Target Name	phospho-Nrf2 (Ser40)
Other Names	phospho-Nrf2 (Ser40)
Accession No.	Uniprot:Q16236Gene ID:4780
Uniprot	Q16236
GeneID	4780
Calculated MW	68
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only