

## PCDH7 Polyclonal Antibody FITC Conjugated

Catalog No: #C91435FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | PCDH7 Polyclonal Antibody FITC Conjugated                                       |
| Host Species       | Rabbit  |
| Clonality          | Polyclonal  |
| Purification       | affinity purified by Protein A  |
| Applications       | Flow-Cyt ICC IF   |
| Species Reactivity | Hu Ms Rat   |
| Immunogen Type     | 31-130/1069   |
| Conjugates         | FITC  |
| Target Name        | PCDH7   |
| Other Names        | PCDH7   |
| Accession No.      | Uniprot:O60245Gene ID:5099  |
| Uniprot            | O60245  |
| GeneID             | 5099  |
| Calculated MW      | 113   |
| Formulation        | Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4        |
| Storage            | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

## Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only