

ALPK3 Polyclonal Antibody FITC Conjugated

Catalog No: #C91436FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ALPK3 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu
Immunogen Type	1821-1907/1907
Conjugates	FITC
Target Name	ALPK3
Other Names	ALPK3
Accession No.	Uniprot:Q96L96Gene ID:57538
Uniprot	Q96L96
GeneID	57538
Calculated MW	201
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only