

# Integrin Alpha V / Beta 1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91440Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Integrin Alpha V / Beta 1 Polyclonal Antibody Biotin Conjugated                 |
| Host Species       | Rabbit  |
| Clonality          | Polyclonal  |
| Purification       | affinity purified by Protein A  |
| Applications       | IHC-P IHC-F IF  |
| Species Reactivity | Hu Ms Rat   |
| Conjugates         | Bio   |
| Target Name        | Integrin Alpha V / Beta 1   |
| Other Names        | Integrin Alpha V / Beta 1   |
| Accession No.      | Uniprot:P06756Gene ID:3685  |
| Uniprot            | P06756  |
| GeneID             | 3685  |
| Calculated MW      | 114   |
| Formulation        | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.                |
| Storage            | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only