

# Integrin Alpha V / Beta 1 Polyclonal Antibody PE-Cy7 Conjugated

Catalog No: #C91440PE-Cy7

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Integrin Alpha V / Beta 1 Polyclonal Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Conjugates	PE-Cy7
Target Name	Integrin Alpha V / Beta 1
Other Names	Integrin Alpha V / Beta 1
Accession No.	Uniprot:P06756Gene ID:3685
Uniprot	P06756
GeneID	3685
Calculated MW	114
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only