

## ALPPL2 Polyclonal Antibody PE Conjugated

Catalog No: #C91442PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ALPPL2 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu
Immunogen Type	301-400/532
Conjugates	PE
Target Name	ALPPL2
Other Names	ALPPL2
Accession No.	Uniprot:P10696Gene ID:251
Uniprot	P10696
GeneID	251
Calculated MW	53
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only