

QPRT/QAPRTase Polyclonal Antibody HRP Conjugated



Catalog No: #C91453HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	QPRT/QAPRTase Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	21-120/297
Conjugates	HRP
Target Name	QPRT/QAPRTase
Other Names	QPRT/QAPRTase
Accession No.	Uniprot:Q15274Gene ID:23475
Uniprot	Q15274
GeneID	23475
Calculated MW	31
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only