

CD31 Polyclonal Antibody PE Conjugated

Catalog No: #C91454PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | CD31 Polyclonal Antibody PE Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | Flow-Cyt ICC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 251-350/727 |
| Conjugates | PE |
| Target Name | CD31 |
| Other Names | CD31 |
| Accession No. | Uniprot:Q08481Gene ID:18613 |
| Uniprot | Q08481 |
| GeneID | 18613 |
| Calculated MW | 78 |
| Formulation | Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4 |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only