

EAAT2 Polyclonal Antibody PE Conjugated

Catalog No: #C91456PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EAAT2 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/574
Conjugates	PE
Target Name	EAAT2
Other Names	EAAT2
Accession No.	Uniprot:P43004Gene ID:6506
Uniprot	P43004
GeneID	6506
Calculated MW	62
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only