

# Myelin-oligodendrocyte glycoprotein Polyclonal Antibody FITC Conjugated

Catalog No: #C91457FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Myelin-oligodendrocyte glycoprotein Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	35-55/247
Conjugates	FITC
Target Name	Myelin-oligodendrocyte glycoprotein
Other Names	Myelin-oligodendrocyte glycoprotein
Accession No.	Uniprot:Q16653Gene ID:4340
Uniprot	Q16653
GeneID	4340
Calculated MW	24
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only