

NFkB p65 Polyclonal Antibody APC Conjugated

Catalog No: #C91459APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NFkB p65 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/551
Conjugates	APC
Target Name	NFkB p65
Other Names	NFkB p65
Accession No.	Uniprot:Q04206Gene ID:5970
Uniprot	Q04206
GeneID	5970
Calculated MW	61
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only