

HIF-1 Alpha Polyclonal Antibody Biotin Conjugated

Catalog No: #C91464Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HIF-1 Alpha Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	661-760/826
Conjugates	Bio
Target Name	HIF-1 Alpha
Other Names	HIF-1 Alpha
Accession No.	Uniprot:Q16665Gene ID:3091
Uniprot	Q16665
GeneID	3091
Calculated MW	92
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only